ÁCIDO ÚRICO ENZIMÁTICO

CÓDIGO	VOLUME (mL)	Nº. TESTES
100/100-100	100	330
100/100-200	200	660

TEST METHODOLOGY	Enzimático
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	510 / 630
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	300
READING TIME (SEC.)	20
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	1
NUMBER OF CUVETTE WASHES	1
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	250
CURVE ACCEPTANCE (%)	20
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	AURI
SAMPLE (µL)	6
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	5.00
TEST LIMIT	20
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	99999
NORMAL VALUES	
MAN	2.5 7.0
WOMAN	1.5 6.0
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	AURI
SAMPLE (µL)	6
PRÉ-DILUTION	10
DILUTION	
FACTOR	5.00
TEST LIMIT	20
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	99999
Normal Values	
Man	250 - 750 mg/24 h
Woman	250 - 750 mg/24 h
Child	-

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TEST METHODOLOGY	VBC
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	630
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	240
READING TIME (SEC.)	10
UNITS -SERUM	g/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	AUB
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	6.0
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	3.5 5.5
WOMAN	1.5 5.5
CHILD	3.8 5.4
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	CNPG
KIND OF PROCESS	Kinetc
TYPE	Linear
FILTERS	405
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	60
READING TIME (SEC.)	240
UNITS -SERUM	U/L
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	AMI
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	7000
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	25 125
WOMAN	25 125
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	DMSO
KIND OF PROCESS	End Point 2 Points
TYPE	Linear
FILTERS	546/630
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	9
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	30 240
READING TIME (SEC.)	20
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	1/1
NUMBER OF CUVETTE WASHES	1
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	300
CURVE ACCEPTANCE (%)	20
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	BILD
SAMPLE (µL)	15
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	10.00
TEST LIMIT	15
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	99999
NORMAL VALUES	
MAN	0 0.25
WOMAN	0 0.25
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	COLOR
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	578
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	30
READING TIME (SEC.)	10
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	CACRE
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	20
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	8.4 10.6
WOMAN	8.4 10.6
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	ARSENAZO
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	630
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	180
READING TIME (SEC.)	10
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	CAARS
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	25
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	8.5 10.5
WOMAN	8.5 10.5
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	TIOCIANATO Hg
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	510
REACTION DIRECTION	Increase
REAGENT #1 (ML)	400
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	180
READING TIME (SEC.)	10
UNITS -SERUM	mEq/L
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1000
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	CLO
SAMPLE (µL)	2
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	150
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	98 100
WOMAN	98 100
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	CLO
SAMPLE (µL)	2
PRÉ-DILUTION	3.00
DILUTION	
FACTOR	2.00
TEST LIMIT	150
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	700
Normal Values	
Man	170 250 mEq/24h
Woman	170 250 mEq/24h
Child	

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TEST METHODOLOGY	IMUNOINIBIÇÃO
KIND OF PROCESS	Kinetc
TYPE	Linear
FILTERS	340
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	300
READING TIME (SEC.)	600
UNITS -SERUM	U/L
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	CKMB
SAMPLE (µL)	12
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	600
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	0-25
WOMAN	0-25
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	CINÉTICO UV
KIND OF PROCESS	Kinetc
TYPE	Linear
FILTERS	340
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	120
READING TIME (SEC.)	300
UNITS -SERUM	U/L
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	CKNAC
SAMPLE (µL)	6
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	2000
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	24 195
WOMAN	24 170
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	ENZIMÁTICO
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	510
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	300
READING TIME (SEC.)	20
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	1/1
NUMBER OF CUVETTE WASHES	1
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	250
CURVE ACCEPTANCE (%)	20
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	COL
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	800
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	99999
NORMAL VALUES	
MAN	0 200
WOMAN	0 200
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	ENZIMÁTICO DIRETO
KIND OF PROCESS	END POINT 2 PT
TYPE	Linear
FILTERS	578
REACTION DIRECTION	Increase
REAGENT #1 (ML)	225
REAGENT #2 (ML)	75
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	300/0
READING TIME (SEC.)	300
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	HDLD
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	200
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	>45
WOMAN	>45
CHILD	>40
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	PRECIPITAÇÃO
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	510/630
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	300
READING TIME (SEC.)	20
UNITS -SERUM	mg/Dl
UNITS -URINE	
NUMBER OF NEEDLE WASHES	1/1
NUMBER OF CUVETTE WASHES	1
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	250
CURVE ACCEPTANCE (%)	20
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	HDLP
SAMPLE (µL)	30
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	1.00
TEST LIMIT	200
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	99999
NORMAL VALUES	
MAN	40 59
WOMAN	40 59
CHILD	>40
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	PICRATO
KIND OF PROCESS	Fixed Time
TYPE	Linear
FILTERS	510
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	60
READING TIME (SEC.)	120
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	CREA
SAMPLE (µL)	30
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	10
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	99999
NORMAL VALUES	
MAN	0.4 1.4
WOMAN	0.4 1.4
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	PICRATO
KIND OF PROCESS	Fixed Time
TYPE	Linear
FILTERS	510
REACTION DIRECTION	Increase
REAGENT #1 (ML)	240
REAGENT #2 (ML)	60
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	60
READING TIME (SEC.)	120
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	CREA
SAMPLE (µL)	30
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	10
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	0.4 1.4
WOMAN	0.4 1.4
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	CROMAZUROL B
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	630
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	300
READING TIME (SEC.)	10
UNITS -SERUM	µg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	FER
SAMPLE (µL)	15
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	400
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	0.4 1.4
WOMAN	0.4 1.4
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	AMP
KIND OF PROCESS	Kinetc
TYPE	Linear
FILTERS	405
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	60
READING TIME (SEC.)	180
UNITS -SERUM	U/L
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	FALC
SAMPLE (µL)	6
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	700
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	0 240
WOMAN	0 270
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	FOSFOMOLIBDATO
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	340
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	0
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	300
READING TIME (SEC.)	10
UNITS -SERUM	mg/dL
UNITS -URINE	mg/dL
NUMBER OF NEEDLE WASHES	2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	FOSF
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	20
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	2000
NORMAL VALUES	
MAN	4.0 7.0
WOMAN	4.0 7.0
CHILD	2.5 5.0
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	FOSF
SAMPLE (µL)	3
PRÉ-DILUTION	20.00
DILUTION	
FACTOR	1.0
TEST LIMIT	20
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	Inactive
Normal Values	
Man	0.0 0.0
Woman	0.0 0.0
Child	0.0 0.0

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TEST METHODOLOGY	NITROBLUETETRAZOLIUM
KIND OF PROCESS	Fixed Time
TYPE	Linear
FILTERS	546
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	600
READING TIME (SEC.)	300
UNITS -SERUM	µmol/L
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	FRUT
SAMPLE (µL)	15
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	800
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	205 285
WOMAN	205 285
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	GLICILGLICINA
KIND OF PROCESS	Kinetc
TYPE	Linear
FILTERS	405
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	60
READING TIME (SEC.)	180
UNITS -SERUM	U/L
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	GGT
SAMPLE (µL)	15
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	250
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	11 50
WOMAN	7 32
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	ENZIMÁTICO
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	510
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	300
READING TIME (SEC.)	10
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	1
NUMBER OF CUVETTE WASHES	1
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	250
CURVE ACCEPTANCE (%)	20
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	GLI
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	500
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	70 99
WOMAN	70 99
CHILD	70 99
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	CINÉTICO UV
KIND OF PROCESS	Kinetc
TYPE	Linear
FILTERS	340
REACTION DIRECTION	Decrease
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	60
READING TIME (SEC.)	210
UNITS -SERUM	U/L
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	GOT
SAMPLE (µL)	30
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	440
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	0-37
WOMAN	0-31
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	ENZIMÁTICO
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	540
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	60
READING TIME (SEC.)	210
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	1
NUMBER OF CUVETTE WASHES	1
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	LAC
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	120
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	4.5 19.8
WOMAN	4.5 19.8
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	PIRUVATO
KIND OF PROCESS	Kinetc
TYPE	Linear
FILTERS	340
REACTION DIRECTION	Decrease
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	30
READING TIME (SEC.)	180
UNITS -SERUM	U/L
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	LDH
SAMPLE (µL)	6
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	2000
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	225 450
WOMAN	225 450
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	AZUL DE XILIDILA
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	510
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	10
INCUBATION TIME (SEC.)	90
READING TIME (SEC.)	20
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1000
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	MG
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	3,5
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	99999
NORMAL VALUES	
MAN	1.9 2.5
WOMAN	1.9 2.5
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	AZUL DE BIURETO
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	540
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	400
READING TIME (SEC.)	10
UNITS -SERUM	g/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	PROT
SAMPLE (µL)	6
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	12
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
NORMAL VALUES	
MAN	6.5 8.0
WOMAN	6.5 8.0
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	PIROGALOL
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	578
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	600
READING TIME (SEC.)	10
UNITS -SERUM	
UNITS -URINE	mg/dL
NUMBER OF NEEDLE WASHES	2/2
NUMBER OF CUVETTE WASHES	2
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
NORMAL VALUES	
MAN	
WOMAN	
CHILD	
CALIBRATOR	
Std. 1	
Std. 2	

URINE	
NAME	PROTUR
SAMPLE (µL)	6
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	2.00
TEST LIMIT	3000
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	500
Normal Values	
Man	10 150
Woman	10 150
Child	10 150

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TEST METHODOLOGY	ENZIMÁTICO
KIND OF PROCESS	End Point
TYPE	Linear
FILTERS	510
REACTION DIRECTION	Increase
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	300
READING TIME (SEC.)	20
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	1
NUMBER OF CUVETTE WASHES	1
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	500
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	TRI
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	5.00
TEST LIMIT	1100
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	99999
NORMAL VALUES	
MAN	60 150
WOMAN	60 150
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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TEST METHODOLOGY	CINÉTICO UV
KIND OF PROCESS	Fixed Time
TYPE	Linear
FILTERS	340
REACTION DIRECTION	Decrease
REAGENT #1 (ML)	300
REAGENT #2 (ML)	-
SAMPLE STARTER	Inactive
DELAY TIME (SEC.)	0
INCUBATION TIME (SEC.)	60
READING TIME (SEC.)	120
UNITS -SERUM	mg/dL
UNITS -URINE	
NUMBER OF NEEDLE WASHES	1/1
NUMBER OF CUVETTE WASHES	1
DYNAMIC BLANK	Inactive
REAGENT BLANK	Every Run
REAGENT LIMIT (mABS)	1000
CURVE ACCEPTANCE (%)	100
INSTRUMENT FACTOR	0.000
SHIFT	0.000

SERUM	
NAME	UREIA
SAMPLE (µL)	3
PRÉ-DILUTION	1.00
DILUTION	
FACTOR	10.00
TEST LIMIT	300
INITIAL ABS	#
FINAL ABS	#
MAX ABS DELTA	1000
NORMAL VALUES	
MAN	15 40
WOMAN	15 40
CHILD	-
CALIBRATOR	
Std. 1	0
Std. 2	*

URINE	
NAME	
SAMPLE (µL)	
PRÉ-DILUTION	
DILUTION	
FACTOR	
TEST LIMIT	
INITIAL ABS	
FINAL ABS	
MAX ABS DELTA	
Normal Values	
Man	
Woman	
Child	

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