ÁCIDO ÚRICO ENZIMÁTICO

CÓDIGO	VOLUME (mL)	Nº. TESTES
100/100-100	100	380
100/100-200	200	760

GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	500 nm
DECIMAL POSITION	2
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	5 µL
DILUTION NAME	H2O
VOLUME	30 µL
REAGENT CYCLE	1
VOLUME	250 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	20.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	T1
LAST	18
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	600 nm
DECIMAL POSITION	2
UNIT	g/dL
ANALYSIS	
POST DIL. FACTOR	NO
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	2 µL
DILUTION NAME	H2O
VOLUME	30 µL
REAGENT CYCLE	1
VOLUME	300 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	OFF
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 g/dL
HIGH	6.0 g/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	CB
LAST	3
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	#
CS2 POS	#
CS3 POS	#

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	FACTOR
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	405 nm
DECIMAL POSITION	1
UNIT	U/L
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	3 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	250 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	NO
HIGH	1000
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	KINETIC
READING FIRST	8
LAST	15
REACTION LIMIT	NO
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
FACTOR	10118
STANDARD POS	*
CONTROL	
CS1 POS	#
CS2 POS	#
CS3 POS	#

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S-SR1
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	550 nm
DECIMAL POSITION	2
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	NO
CONC. FACTOR	NO
SAMPLE CYCLE	2
VOLUME	25 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	375 µL
START R1 CYCLE	1
VOLUME	15 µL
DILUENT NAME	H2O
VOLUME	10 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	15.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	1
LAST	10
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	#
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	600 nm
DECIMAL POSITION	1
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	NO
CONC. FACTOR	NO
SAMPLE CYCLE	2
VOLUME	3 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	300 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	NO
HIGH	20.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	1
LAST	6
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	#
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	FACTOR
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	340 nm
DECIMAL POSITION	1
UNIT	U/L
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	4 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	190 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 U/L
HIGH	2000.0 U/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	KINETIC
READING FIRST	5
LAST	13
REACTION LIMIT	0.2500 ΔA
POINT	1
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
FACTOR	17976
STANDARD POS	*
CONTROL	
CS1 POS	#
CS2 POS	#
CS3 POS	#

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	FACTOR
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	340 nm
DECIMAL POSITION	1
UNIT	U/L
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	8 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	190 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	NO
HIGH	600.0 U/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	KINETIC
READING FIRST	5
LAST	18
REACTION LIMIT	0.3500 ΔA
POINT	1
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
FACTOR	4676
STANDARD POS	*
CONTROL	
CS1 POS	#
CS2 POS	#
CS3 POS	#

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	500 nm
DECIMAL POSITION	1
UNIT	mEq/L
ANALYSIS	
POST DIL. FACTOR	2.00
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	2 µL
DILUTION NAME	H2O
VOLUME	20 µL
REAGENT CYCLE	1
VOLUME	400 µL
CALCULATION	
SAMPLE LIMIT	0.1000 ΔA
POINT	T1
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	70.0 mEq/L
HIGH	150.0 mEq/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	CB
LAST	10
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	500 nm
DECIMAL POSITION	1
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	2.00
CONC. FACTOR	NO
SAMPLE CYCLE	2
VOLUME	3 µL
DILUTION NAME	H2O
VOLUME	30 µL
REAGENT CYCLE	1
VOLUME	300 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	800.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	1
LAST	12
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	500 nm
DECIMAL POSITION	1
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	2.00
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	20 µL
DILUTION NAME	H2O
VOLUME	30 µL
REAGENT CYCLE	1
VOLUME	200 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	150.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	CB
LAST	18
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S-RS1
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	600 nm
DECIMAL POSITION	1
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	NO
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	2.5 µL
DILUTION NAME	H2O
VOLUME	20 µL
REAGENT CYCLE	1
VOLUME	190 µL
START R1 CICLE	12
VOLUME	65
DILUTION NAME	H2O
VOLUME	10
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	NO
HIGH	180.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	11
LAST	24
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	CALIBRATOR
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	500 nm
DECIMAL POSITION	2
UNIT	mg/dL
ANALYSIS	
SAMPLE DIL NAME	H2O
POST DIL. FACTOR	2.00
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	20 µL
DILUTION NAME	H2O
VOLUME	30 ML
REAGENT CYCLE	1
VOLUME	200 µL
CALCULATION	
SAMPLE LIMIT	0.3600 ΔA
POINT	T1
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	-0.2500
TEST RANGE LOW	0.20 mg/dL
HIGH	10.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	2
LAST	9
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	-0.0090 ΔA
HIGH	-0.2800 ΔA
REAGENT RANGE LOW	-0.0090 ΔA
HIGH	-0.0200 ΔA
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S-SR1
CALIBRATION MODE	CALIBRATOR
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	500 nm
DECIMAL POSITION	2
UNIT	mg/dL
ANALYSIS	
SAMPLE DIL NAME	H2O
POST DIL. FACTOR	2.00
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	20 µL
DILUTION NAME	H2O
VOLUME	30 µL
REAGENT CYCLE	1
VOLUME	160 µL
START R1 CYCLE	2
VOL	40 µL
DIL	0.0 µL
CALCULATION	
SAMPLE LIMIT	0.3600 ΔA
POINT	T1
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	-0.25000
TEST RANGE LOW	0.20 mg/dL
HIGH	10.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	3
LAST	10
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	-0.0090 ΔA
HIGH	-0.2800 ΔA
REAGENT RANGE LOW	-0.0090 ΔA
HIGH	-0.0200 ΔA
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	CALIBRATOR
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	600 nm
DECIMAL POSITION	1
UNIT	µg/dL
ANALYSIS	
POST DIL. FACTOR	NO
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	10 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	200 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 µg/dL
HIGH	500.0 µg/dL
NORM RANGE LOW	37 µg/dL
HIGH	148 µg/dL
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	3
LAST	27
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	5.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	CALIBRATOR
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	550 nm
DECIMAL POSITION	1
UNIT	µg/dL
ANALYSIS	
POST DIL. FACTOR	NO
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	50 µL
DILUTION NAME	H2O
VOLUME	30 µL
REAGENT CYCLE	1
VOLUME	160 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 µg/dL
HIGH	400.0 µg/Dl
NORM RANGE LOW	45 µg/dL
HIGH	150 µg/dL
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	3
LAST	27
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	5.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	FACTOR
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	405 nm
DECIMAL POSITION	1
UNIT	U/L
ANALYSIS	
POST DIL. FACTOR	5
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	5 µL
DILUTION NAME	H2O
VOLUME	25 µL
REAGENT CYCLE	1
VOLUME	225 µL
CALCULATION	
SAMPLE LIMIT	0.2000 ΔA
POINT	T1
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 U/L
HIGH	700 U/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	KINETIC
READING FIRST	4
LAST	12
REACTION LIMIT	0.2500 ΔA
POINT	T1
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
FACTOR	4594
STANDARD POS	*
CONTROL	
CS1 POS	#
CS2 POS	#
CS3 POS	#

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	340 nm
DECIMAL POSITION	2
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	2
VOLUME	4.0 µL
DILUTION NAME	H2O
VOLUME	20 µL
REAGENT CYCLE	1
VOLUME	250 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	15.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	1
LAST	6
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOP AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	550 nm
DECIMAL POSITION	2
UNIT	µmol/L
ANALYSIS	
POST DIL. FACTOR	NO
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	15 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	300 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 µmol/L
HIGH	800.0 µmol/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	25
LAST	37
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	FACTOR
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	405 nm
DECIMAL POSITION	1
UNIT	U/L
ANALYSIS	
POST DIL. FACTOR	5
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	25 µL
DILUTION NAME	H2O
VOLUME	25 µL
REAGENT CYCLE	1
VOLUME	225 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 U/L
HIGH	250 U/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	KINETIC
READING FIRST	4
LAST	8
REACTION LIMIT	0.4300 ΔA
POINT	T1
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
FACTOR	2655
STANDARD POS	*
CONTROL	
CS1 POS	#
CS2 POS	#
CS3 POS	#

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	500 nm
DECIMAL POSITION	1
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	2.5 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	250 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	500.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	1
LAST	20
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	NO
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	FACTOR
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	340 nm
DECIMAL POSITION	1
UNIT	U/L
ANALYSIS	
POST DIL. FACTOR	8
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	15 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	150 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	DECREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	NO
HIGH	440 U/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	KINSEARCH
READING FIRST	3
LAST	11
REACTION LIMIT	NO
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
FACTOR	4286
STANDARD POS	*
CONTROL	
CS1 POS	#
CS2 POS	#
CS3 POS	#

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	FACTOR
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	340 nm
DECIMAL POSITION	1
UNIT	U/L
ANALYSIS	
POST DIL. FACTOR	8
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	15 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	150 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	DECREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	NO
HIGH	350 U/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	KINSEARCH
READING FIRST	3
LAST	11
REACTION LIMIT	NO
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
FACTOR	4477
STANDARD POS	*
CONTROL	
CS1 POS	#
CS2 POS	#
CS3 POS	#

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	550 nm
DECIMAL POSITION	1
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	3 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	300 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	120.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	1
LAST	12
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	FACTOR
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	340 nm
DECIMAL POSITION	1
UNIT	U/L
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	5 µL
DILUTION NAME	H2O
VOLUME	20 µL
REAGENT CYCLE	1
VOLUME	225 µL
CALCULATION	
SAMPLE LIMIT	0.2000 ΔA
POINT	T1
REACTION DIRECTION	DECREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 U/L
HIGH	2000.0 U/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	KINETIC
READING FIRST	3
LAST	11
REACTION LIMIT	0.2000 ΔA
POINT	T1
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
FACTOR	13227
STANDARD POS	*
CONTROL	
CS1 POS	#
CS2 POS	#
CS3 POS	#

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S-RS1
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	550 nm
DECIMAL POSITION	1
UNIT	U/L
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	3 µL
DILUTION NAME	H2O
VOLUME	20 µL
REAGENT CYCLE	1
VOLUME	240 µL
START R1 CYCLE	2
VOLUME	60 µL
DILUENT NAME	H2O
VOLUME	0.0 ml
CALCULATION	
SAMPLE LIMIT	0.3600 ΔA
POINT	T1
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 U/L
HIGH	00.0 U/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	6
LAST	12
REACTION	NO
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	500 nm
DECIMAL POSITION	2
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	23 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	300 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	3.5 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	CB
LAST	6
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	578 nm
DECIMAL POSITION	2
UNIT	mmol/L
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	2
VOLUME	30 µL
DILUTION NAME	H2O
VOLUME	30 µL
REAGENT CYCLE	1
VOLUME	300 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mmol/L
HIGH	10.0 mmol/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	T1
LAST	8
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	550 nm
DECIMAL POSITION	2
UNIT	g/dL
ANALYSIS	
POST DIL. FACTOR	2
CONC. FACTOR	NO
SAMPLE CYCLE	2
VOLUME	3 µL
DILUTION NAME	H2O
VOLUME	30 µL
REAGENT CYCLE	1
VOLUME	250 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 g/dL
HIGH	12.0 g/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	T1
LAST	T8
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	NO BLANK
CLEANER	NO
WAVELENGTH	600 nm
DECIMAL POSITION	1
UNIT	mg/L
ANALYSIS	
POST DIL. FACTOR	NO
CONC. FACTOR	3
SAMPLE CYCLE	2
VOLUME	5 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	300 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	30.0 mg/L
HIGH	3000.0 mg/L
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	1
LAST	12
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	BEFORE
WAVELENGTH	500 nm
DECIMAL POSITION	1
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	3
CONC. FACTOR	NO
SAMPLE CYCLE	2
VOLUME	3 µL
DILUTION NAME	H2O
VOLUME	10 µL
REAGENT CYCLE	1
VOLUME	250 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	INCREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	800.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	1
LAST	18
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
TIME	NO
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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GENERAL	
MEASUREMENT MODE	ABSORB
REACTION MODE	R-S
CALIBRATION MODE	SLOPE AVG
REAGENT BLANK	REAG/DIL
CLEANER	NO
WAVELENGTH	340 nm
DECIMAL POSITION	1
UNIT	mg/dL
ANALYSIS	
POST DIL. FACTOR	3
CONC. FACTOR	NO
SAMPLE CYCLE	1
VOLUME	3 µL
DILUTION NAME	H2O
VOLUME	30 µL
REAGENT CYCLE	1
VOLUME	220 µL
CALCULATION	
SAMPLE LIMIT	NO
REACTION DIRECTION	DECREASE
CHECK	ON
CONVERSION FACTOR	1.000
OFFSET	0.000
TEST RANGE LOW	0.0 mg/dL
HIGH	250.0 mg/dL
NORM RANGE LOW	NO
HIGH	NO
NUMBER OF STEPS	1
CALCULATION STEPS A	ENDPOINT
READING FIRST	2
LAST	8
CALIBRATION	
CALIB. INTERVAL	ON REQUEST
BLANK RANGE LOW	NO
HIGH	NO
REAGENT RANGE LOW	NO
HIGH	NO
STANDARD POS	#
STD-1	*
STD-2	NO
STD-3	NO
REPLICATE	DUPL
DEVIATION	10.0 %
CONTROL	
CS1 POS	*
CS2 POS	*
CS3 POS	*

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