ÁCIDO ÚRICO ENZIMÁTICO

CÓDIGO	VOLUME (mL)	Nº. TESTES
100/100-100	100	500
100/100-200	200	1000

TEST NAME	AUR		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	4	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	520		SEC				
METHOD		END						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	20.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	mg/dL		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	ALB	TYPE	SERUM	OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl
REAGENTS	R1 VOLUME	400	µl	DILUTION	0	µl
	R2 VOLUME		µl	DILUTION	0	µl
WAVELENGTH	PRI	660		SEC		
METHOD		END				
REACTION SLOP		+				
MEASURING POIN 1	FIRST	0		LAST	10	
MEASURING POIN 2	FIRST			LAST		
LINEARITY		%				
N0-LAG-TIME						
PRE-DILUTION RATE						
	MIN OD			MAX OD		
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50	
	LAST L	-0.100		LAST H	2.50	
DYNAMIC RANGE	L	1.0		H	6.0	
CORRELATION FACTOR	A	1		B	0	
ON-BOARD STABILITY PERIOD		30 DAYS				
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#
			AGE L		AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH	
	1	#	#	#	#	
	2	#	#	#	#	
	3	#	#	#	#	
	4	#	#	#	#	
	5	#	#	#	#	
	6	#	#	#	#	
	7	#	#	#	#	
PANIC VALUE	L	#	H	#	UNIT	g/dL
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H	
POINT 1	*		*			
POINT 2	*		*			
POINT 3	*		*			
POINT 4	*		*			
POINT 5	*		*			
POINT 6	*		*			
POINT 7	*		*			
1-POINT CAL. POINT						
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD		

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TEST NAME	AMIL		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	410		SEC				
METHOD		RATE						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	10			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	1.0		H	2000			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	U/L		
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	AIGP		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD		END						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	1.0		H	150			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#	AGE L	AGE H
			SEX	YEAR	MONTH	YEAR	MONTH	
			1	#	#	#	#	#
			2	#	#	#	#	#
			3	#	#	#	#	#
			4	#	#	#	#	#
			5	#	#	#	#	#
			6	#	#	#	#	#
			7	#	#	#	#	#
PANIC VALUE	L	#	H	#	UNIT	mg/dL		
CALIBRATION TYPE	AB		FORMULA	Y=AX+B	COUNTS	#		
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD				

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TEST NAME	BILD		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	10	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME	6	µl	DILUTION	0	µl		
WAVELENGTH	PRI	540		SEC				
METHOD		END						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	15.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	mg/dL		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B		COUNTS	#		
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR			CALIBRATION STABILITY PERIOD					

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TEST NAME	BILT		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	10	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME	6	µl	DILUTION	0	µl		
WAVELENGTH	PRI	540		SEC				
METHOD		END						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	15.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#	AGE L	AGE H
			SEX	YEAR	MONTH	YEAR	MONTH	
			1	#	#	#	#	#
			2	#	#	#	#	#
			3	#	#	#	#	#
			4	#	#	#	#	#
			5	#	#	#	#	#
			6	#	#	#	#	#
			7	#	#	#	#	#
PANIC VALUE	L	#	H	#	UNIT	mg/dL		
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR	FACTOR			
				OD-L	OD-H			
	POINT 1	*		*				
	POINT 2	*		*				
	POINT 3	*		*				
	POINT 4	*		*				
	POINT 5	*		*				
	POINT 6	*		*				
	POINT 7	*		*				
1-POINT CAL. POINT								
MB TYPE FACTOR								
	CALIBRATION STABILITY PERIOD							

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TEST NAME	ASO TURBI		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	540		SEC				
METHOD	FIXED							
REACTION SLOP	+							
MEASURING POIN 1	FIRST	1		LAST	8			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	800.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#	AGE L	AGE H
	SEX	YEAR	MONTH	YEAR	MONTH			
	1	#	#	#	#	#	#	#
	2	#	#	#	#	#	#	#
	3	#	#	#	#	#	#	#
	4	#	#	#	#	#	#	#
	5	#	#	#	#	#	#	#
	6	#	#	#	#	#	#	#
	7	#	#	#	#	#	#	#
PANIC VALUE	L	#		H	#	UNIT	UI/mL	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR			CALIBRATION STABILITY PERIOD					

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TEST NAME	CAL	TYPE	SERUM	OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl
	R2 VOLUME		µl	DILUTION	0	µl
WAVELENGTH	PRI	660		SEC		
METHOD		END				
REACTION SLOP		+				
MEASURING POIN 1	FIRST	0		LAST	27	
MEASURING POIN 2	FIRST			LAST		
LINEARITY		%				
N0-LAG-TIME						
PRE-DILUTION RATE						
	MIN OD			MAX OD		
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50	
	LAST L	-0.100		LAST H	2.50	
DYNAMIC RANGE	L	0.0		H	20.0	
CORRELATION FACTOR	A	1		B	0	
ON-BOARD STABILITY PERIOD		30 DAYS				
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#
			AGE L		AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH	
	1	#	#	#	#	
	2	#	#	#	#	
	3	#	#	#	#	
	4	#	#	#	#	
	5	#	#	#	#	
	6	#	#	#	#	
	7	#	#	#	#	
PANIC VALUE	L	#	H	#	UNIT	mg/dL
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H	
POINT 1	*		*			
POINT 2	*		*			
POINT 3	*		*			
POINT 4	*		*			
POINT 5	*		*			
POINT 6	*		*			
POINT 7	*		*			
1-POINT CAL. POINT						
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD		

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TEST NAME	<input type="text" value="CLO"/>	TYPE	<input type="text" value="SERUM"/>	OPERATION	<input type="text" value="YES"/>	
SAMPLE	VOLUME	<input type="text" value="2"/>	µl	DILUTION	<input type="text" value="0"/>	µl
REAGENTS	R1 VOLUME	<input type="text" value="400"/>	µl	DILUTION	<input type="text" value="0"/>	µl
	R2 VOLUME	<input type="text"/>	µl	DILUTION	<input type="text" value="0"/>	µl
WAVELENGTH	PRI	<input type="text" value="520"/>		SEC	<input type="text"/>	
METHOD		<input type="text" value="END"/>				
REACTION SLOP		<input type="text" value="+"/>				
MEASURING POIN 1	FIRST	<input type="text" value="0"/>		LAST	<input type="text" value="10"/>	
MEASURING POIN 2	FIRST	<input type="text"/>		LAST	<input type="text"/>	
LINEARITY		<input type="text"/>	%			
N0-LAG-TIME		<input type="text"/>				
PRE-DILUTION RATE		<input type="text"/>				
	MIN OD	<input type="text"/>		MAX OD	<input type="text"/>	
REAGENT OD LIMIT	FIRST L	<input type="text" value="-0.100"/>		FIRST H	<input type="text" value="2.50"/>	
	LAST L	<input type="text" value="-0.100"/>		LAST H	<input type="text" value="2.50"/>	
DYNAMIC RANGE		L	<input type="text" value="1.0"/>	H	<input type="text" value="150.0"/>	
CORRELATION FACTOR		A	<input type="text" value="1"/>	B	<input type="text" value="0"/>	
ON-BOARD STABILITY PERIOD		<input type="text" value="30 DAYS"/>				
NORMAL RANGES	VALUE/FLAG	<input type="text" value="#"/>	LEVE L	<input type="text" value="#"/>	LEVEL H	<input type="text" value="#"/>
			AGE L		AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH	
	1	#	#	#	#	
	2	#	#	#	#	
	3	#	#	#	#	
	4	#	#	#	#	
	5	#	#	#	#	
	6	#	#	#	#	
	7	#	#	#	#	
PANIC VALUE	L	<input type="text" value="#"/>	H	<input type="text" value="#"/>	UNIT	<input type="text" value="mEq/L"/>
CALIBRATION TYPE	<input type="text" value="AB"/>	FORMULA	<input type="text" value="Y=AX+B"/>	COUNTS	<input type="text" value="#"/>	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H	
POINT 1	*		*			
POINT 2	*		*			
POINT 3	*		*			
POINT 4	*		*			
POINT 5	*		*			
POINT 6	*		*			
POINT 7	*		*			
1-POINT CAL. POINT	<input type="text"/>					
MB TYPE FACTOR	<input type="text"/>			CALIBRATION STABILITY PERIOD	<input type="text"/>	

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TEST NAME	CK MB		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	8	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD	RATE							
REACTION SLOP	+							
MEASURING POIN 1	FIRST	22		LAST	34			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	600.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#		H	#	UNIT	U/L	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR								
						CALIBRATION STABILITY PERIOD		

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TEST NAME	CK MB		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	8	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD	RATE							
REACTION SLOP		+						
MEASURING POIN 1	FIRST	12		LAST	24			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	600.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	U/L		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	CK T		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	4	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD	RATE							
REACTION SLOP		+						
MEASURING POIN 1	FIRST	17		LAST	29			
MEASURING POIN 2	FIRST							
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	2000.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#	AGE L	AGE H
	SEX	YEAR	MONTH	YEAR	MONTH			
	1	#	#	#	#	#	#	#
	2	#	#	#	#	#	#	#
	3	#	#	#	#	#	#	#
	4	#	#	#	#	#	#	#
	5	#	#	#	#	#	#	#
	6	#	#	#	#	#	#	#
	7	#	#	#	#	#	#	#
PANIC VALUE	L	#		H	#	UNIT	U/L	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR								
	CALIBRATION STABILITY PERIOD							

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TEST NAME	CK T		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	4	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD	RATE							
REACTION SLOP		+						
MEASURING POIN 1	FIRST	12		LAST	24			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	2000.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	U/L		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	COLT		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	520		SEC				
METHOD		END						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	800.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	mg/dL		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	HDL D		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	3	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	225	µl	DILUTION	0	µl		
	R2 VOLUME	75	µl	DILUTION	0	µl		
WAVELENGTH	PRI	570		SEC				
METHOD	END							
REACTION SLOP	+							
MEASURING POIN 1	FIRST	10		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	150.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#		H	#	UNIT	mg/dL	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR	FACTOR			
POINT 1	*		*	OD-L	OD-H			
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR			CALIBRATION STABILITY PERIOD					

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TEST NAME	CREA	TYPE	SERUM	OPERATION	YES	
SAMPLE	VOLUME	20	µl	DILUTION	0	µl
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl
	R2 VOLUME		µl	DILUTION	0	µl
WAVELENGTH	PRI	520		SEC		
METHOD		FIXED				
REACTION SLOP		+				
MEASURING POIN 1	FIRST	4		LAST	12	
MEASURING POIN 2	FIRST			LAST		
LINEARITY		%				
N0-LAG-TIME						
PRE-DILUTION RATE						
	MIN OD			MAX OD		
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50	
	LAST L	-0.100		LAST H	2.50	
DYNAMIC RANGE	L	0.0		H	10.0	
CORRELATION FACTOR	A	1		B	0	
ON-BOARD STABILITY PERIOD		30 DAYS				
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#
			AGE L		AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH	
	1	#	#	#	#	
	2	#	#	#	#	
	3	#	#	#	#	
	4	#	#	#	#	
	5	#	#	#	#	
	6	#	#	#	#	
	7	#	#	#	#	
PANIC VALUE	L	#	H	#	UNIT	mg/dL
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H	
POINT 1	*		*			
POINT 2	*		*			
POINT 3	*		*			
POINT 4	*		*			
POINT 5	*		*			
POINT 6	*		*			
POINT 7	*		*			
1-POINT CAL. POINT						
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD		

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TEST NAME	CREA	TYPE	SERUM	OPERATION	YES	
SAMPLE	VOLUME	20	µl	DILUTION	0	µl
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl
	R2 VOLUME	40	µl	DILUTION	0	µl
WAVELENGTH	PRI	520		SEC		
METHOD		FIXED				
REACTION SLOP		+				
MEASURING POIN 1	FIRST	10		LAST	18	
MEASURING POIN 2	FIRST			LAST		
LINEARITY		%				
N0-LAG-TIME						
PRE-DILUTION RATE						
	MIN OD			MAX OD		
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50	
	LAST L	-0.100		LAST H	2.50	
DYNAMIC RANGE	L	0.0		H	10.0	
CORRELATION FACTOR	A	1		B	0	
ON-BOARD STABILITY PERIOD		30 DAYS				
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#
			AGE L		AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH	
	1	#	#	#	#	
	2	#	#	#	#	
	3	#	#	#	#	
	4	#	#	#	#	
	5	#	#	#	#	
	6	#	#	#	#	
	7	#	#	#	#	
PANIC VALUE	L	#	H	#	UNIT	mg/dL
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H	
POINT 1	*		*			
POINT 2	*		*			
POINT 3	*		*			
POINT 4	*		*			
POINT 5	*		*			
POINT 6	*		*			
POINT 7	*		*			
1-POINT CAL. POINT						
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD		

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TEST NAME	FERRITINA		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	18	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	540		SEC				
METHOD	FIXED							
REACTION SLOP	+							
MEASURING POIN 1	FIRST	4		LAST	24			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	600.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#		H	#	UNIT	µg/L	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR	FACTOR			
POINT 1	*		*	OD-L	OD-H			
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR			CALIBRATION STABILITY PERIOD					

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TEST NAME	FE CROM		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	10	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	620		SEC				
METHOD		END						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	500.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	µg/dL		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	FALC		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	4	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	405		SEC				
METHOD		RATE						
REACTION SLOP		+						
MEASURING POIN 1		FIRST	13	LAST	25			
MEASURING POIN 2		FIRST		LAST				
LINEARITY			%					
N0-LAG-TIME								
PRE-DILUTION RATE								
		MIN OD		MAX OD				
REAGENT OD LIMIT		FIRST L	-0.100	FIRST H	2.50			
		LAST L	-0.100	LAST H	2.50			
DYNAMIC RANGE		L	0.0	H	700.0			
CORRELATION FACTOR		A	1	B	0			
ON-BOARD STABILITY PERIOD		30 DAYS						
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
		SEX	YEAR	MONTH	YEAR	MONTH		
	1	#	#	#	#	#		
	2	#	#	#	#	#		
	3	#	#	#	#	#		
	4	#	#	#	#	#		
	5	#	#	#	#	#		
	6	#	#	#	#	#		
	7	#	#	#	#	#		
PANIC VALUE		L	#	H	#	UNIT	U/L	
CALIBRATION TYPE		AB	FORMULA	Y=AX+B	COUNTS	#		
		CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H		
	POINT 1	*		*				
	POINT 2	*		*				
	POINT 3	*		*				
	POINT 4	*		*				
	POINT 5	*		*				
	POINT 6	*		*				
	POINT 7	*		*				
1-POINT CAL. POINT								
MB TYPE FACTOR					CALIBRATION STABILITY PERIOD			

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TEST NAME	FALC		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	4	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	405		SEC				
METHOD		RATE						
REACTION SLOP		+						
MEASURING POIN 1		FIRST	8		LAST	20		
MEASURING POIN 2		FIRST			LAST			
LINEARITY			%					
N0-LAG-TIME								
PRE-DILUTION RATE								
		MIN OD			MAX OD			
REAGENT OD LIMIT		FIRST L	-0.100		FIRST H	2.50		
		LAST L	-0.100		LAST H	2.50		
DYNAMIC RANGE		L	0.0		H	700.0		
CORRELATION FACTOR		A	1		B	0		
ON-BOARD STABILITY PERIOD		30 DAYS						
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
		SEX	YEAR	MONTH	YEAR	MONTH		
	1	#	#	#	#	#		
	2	#	#	#	#	#		
	3	#	#	#	#	#		
	4	#	#	#	#	#		
	5	#	#	#	#	#		
	6	#	#	#	#	#		
	7	#	#	#	#	#		
PANIC VALUE		L	#		H	#	UNIT	U/L
CALIBRATION TYPE		AB		FORMULA	Y=AX+B		COUNTS	#
		CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H		
	POINT 1	*		*				
	POINT 2	*		*				
	POINT 3	*		*				
	POINT 4	*		*				
	POINT 5	*		*				
	POINT 6	*		*				
	POINT 7	*		*				
1-POINT CAL. POINT								
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD				

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TEST NAME	FOSF		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD		END						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	8			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	30.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	mg/dL		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	FRUTO		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	10	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	520		SEC				
METHOD		FIXED						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	20		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	7.00			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	800.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD		30 DAYS						
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#	AGE L	AGE H
			SEX	YEAR	MONTH	YEAR	MONTH	
			1	#	#	#	#	#
			2	#	#	#	#	#
			3	#	#	#	#	#
			4	#	#	#	#	#
			5	#	#	#	#	#
			6	#	#	#	#	#
			7	#	#	#	#	#
PANIC VALUE	L	#	H	#	UNIT	µmol/L		
CALIBRATION TYPE	AB		FORMULA	Y=AX+B	COUNTS	#		
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD				

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TEST NAME	FR TURBI		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	2.1	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	240	µl	DILUTION	0	µl		
	R2 VOLUME	60	µl	DILUTION	0	µl		
WAVELENGTH	PRI	640		SEC				
METHOD		FIXED						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	4		LAST	12			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	160.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	UI/mL		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	GGT		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	20	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	405		SEC				
METHOD		RATE						
REACTION SLOP		+						
MEASURING POIN 1		FIRST	13	LAST	25			
MEASURING POIN 2		FIRST		LAST				
LINEARITY			%					
N0-LAG-TIME								
PRE-DILUTION RATE								
		MIN OD		MAX OD				
REAGENT OD LIMIT		FIRST L	-0.100	FIRST H	2.50			
		LAST L	-0.100	LAST H	2.50			
DYNAMIC RANGE		L	0.0	H	250.0			
CORRELATION FACTOR		A	1	B	0			
ON-BOARD STABILITY PERIOD		30 DAYS						
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#	AGE L	AGE H
			SEX	YEAR	MONTH	YEAR	MONTH	
			1	#	#	#	#	#
			2	#	#	#	#	#
			3	#	#	#	#	#
			4	#	#	#	#	#
			5	#	#	#	#	#
			6	#	#	#	#	#
			7	#	#	#	#	#
PANIC VALUE		L	#	H	#	UNIT	U/L	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR	FACTOR			
	POINT 1	*		*				
	POINT 2	*		*				
	POINT 3	*		*				
	POINT 4	*		*				
	POINT 5	*		*				
	POINT 6	*		*				
	POINT 7	*		*				
1-POINT CAL. POINT								
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD				

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TEST NAME	GGT		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	20	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	405		SEC				
METHOD		RATE						
REACTION SLOP		+						
MEASURING POIN 1		FIRST	8	LAST	20			
MEASURING POIN 2		FIRST		LAST				
LINEARITY			%					
N0-LAG-TIME								
PRE-DILUTION RATE								
		MIN OD		MAX OD				
REAGENT OD LIMIT		FIRST L	-0.100	FIRST H	2.50			
		LAST L	-0.100	LAST H	2.50			
DYNAMIC RANGE		L	0.0	H	250.0			
CORRELATION FACTOR		A	1	B	0			
ON-BOARD STABILITY PERIOD		30 DAYS						
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#	AGE L	AGE H
			SEX	YEAR	MONTH	YEAR	MONTH	
			1	#	#	#	#	#
			2	#	#	#	#	#
			3	#	#	#	#	#
			4	#	#	#	#	#
			5	#	#	#	#	#
			6	#	#	#	#	#
			7	#	#	#	#	#
PANIC VALUE		L	#	H	#	UNIT	U/L	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR	FACTOR			
	POINT 1	*		OD-L	OD-H			
	POINT 2	*						
	POINT 3	*						
	POINT 4	*						
	POINT 5	*						
	POINT 6	*						
	POINT 7	*						
1-POINT CAL. POINT								
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD				

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TEST NAME	GLIC		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	<input type="text" value="2"/>	µl	DILUTION	<input type="text" value="0"/>	µl		
REAGENTS	R1 VOLUME	<input type="text" value="200"/>	µl	DILUTION	<input type="text" value="0"/>	µl		
	R2 VOLUME	<input type="text"/>	µl	DILUTION	<input type="text" value="0"/>	µl		
WAVELENGTH	PRI	<input type="text" value="520"/>		SEC	<input type="text"/>			
METHOD		<input type="text" value="END"/>						
REACTION SLOP		<input type="text" value="+"/>						
MEASURING POIN 1	FIRST	<input type="text" value="0"/>		LAST	<input type="text" value="27"/>			
MEASURING POIN 2	FIRST	<input type="text"/>		LAST	<input type="text"/>			
LINEARITY		<input type="text"/>	%					
N0-LAG-TIME		<input type="text"/>						
PRE-DILUTION RATE		<input type="text"/>						
	MIN OD	<input type="text"/>		MAX OD	<input type="text"/>			
REAGENT OD LIMIT	FIRST L	<input type="text" value="-0.100"/>		FIRST H	<input type="text" value="2.50"/>			
	LAST L	<input type="text" value="-0.100"/>		LAST H	<input type="text" value="2.50"/>			
DYNAMIC RANGE	L	<input type="text" value="0.0"/>		H	<input type="text" value="500.0"/>			
CORRELATION FACTOR	A	<input type="text" value="1"/>		B	<input type="text" value="0"/>			
ON-BOARD STABILITY PERIOD	<input type="text" value="30 DAYS"/>							
NORMAL RANGES	VALUE/FLAG	<input type="text" value="#"/>	LEVE L	<input type="text" value="#"/>	LEVEL H	<input type="text" value="#"/>		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>			
2	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>			
3	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>			
4	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>			
5	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>			
6	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>			
7	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>	<input type="text" value="#"/>			
PANIC VALUE	L	<input type="text" value="#"/>	H	<input type="text" value="#"/>	UNIT	<input type="text" value="mg/dL"/>		
CALIBRATION TYPE	<input type="text" value="AB"/>		FORMULA	<input type="text" value="Y=AX+B"/>		COUNTS	<input type="text" value="#"/>	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>			
POINT 2	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>			
POINT 3	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>			
POINT 4	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>			
POINT 5	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>			
POINT 6	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>			
POINT 7	<input type="text" value="*"/>	<input type="text"/>	<input type="text" value="*"/>	<input type="text"/>	<input type="text"/>			
1-POINT CAL. POINT	<input type="text"/>							
MB TYPE FACTOR	<input type="text"/>							
	CALIBRATION STABILITY PERIOD				<input type="text"/>			

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TEST NAME	TGO		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	20	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD	RATE							
REACTION SLOP		-						
MEASURING POIN 1	FIRST	8		LAST	20			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	440.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#		LEVE L	#		LEVEL H	#
				AGE L			AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH			
	1	#	#	#	#			
	2	#	#	#	#			
	3	#	#	#	#			
	4	#	#	#	#			
	5	#	#	#	#			
	6	#	#	#	#			
	7	#	#	#	#			
PANIC VALUE	L	#		H	#		UNIT	U/L
CALIBRATION TYPE	AB			FORMULA	Y=AX+B		COUNTS	#
	CAL N°	OD	CONC	FACTOR	FACTOR			
	POINT 1	*		OD-L	OD-H			
	POINT 2	*						
	POINT 3	*						
	POINT 4	*						
	POINT 5	*						
	POINT 6	*						
	POINT 7	*						
1-POINT CAL. POINT								
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD				

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TEST NAME	TGO		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	20	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD		RATE						
REACTION SLOP		-						
MEASURING POIN 1		FIRST	13	LAST	25			
MEASURING POIN 2		FIRST		LAST				
LINEARITY			%					
N0-LAG-TIME								
PRE-DILUTION RATE								
		MIN OD		MAX OD				
REAGENT OD LIMIT		FIRST L	-0.100	FIRST H	2.50			
		LAST L	-0.100	LAST H	2.50			
DYNAMIC RANGE		L	0.0	H	440.0			
CORRELATION FACTOR		A	1	B	0			
ON-BOARD STABILITY PERIOD		30 DAYS						
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#	AGE L	AGE H
			SEX	YEAR	MONTH	YEAR	MONTH	
			1	#	#	#	#	#
			2	#	#	#	#	#
			3	#	#	#	#	#
			4	#	#	#	#	#
			5	#	#	#	#	#
			6	#	#	#	#	#
			7	#	#	#	#	#
PANIC VALUE		L	#	H	#	UNIT	U/L	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL Nº	OD	CONC	FACTOR	FACTOR			
	POINT 1	*		OD-L	OD-H			
	POINT 2	*		*				
	POINT 3	*		*				
	POINT 4	*		*				
	POINT 5	*		*				
	POINT 6	*		*				
	POINT 7	*		*				
1-POINT CAL. POINT								
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD				

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TEST NAME	TGP		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	20	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD	RATE							
REACTION SLOP		-						
MEASURING POIN 1	FIRST	8		LAST	20			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	350.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	U/L		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	TGP		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	20	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD		RATE						
REACTION SLOP		-						
MEASURING POIN 1		FIRST	13	LAST	25			
MEASURING POIN 2		FIRST		LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
		MIN OD		MAX OD				
REAGENT OD LIMIT		FIRST L	-0.100	FIRST H	2.50			
		LAST L	-0.100	LAST H	2.50			
DYNAMIC RANGE		L	0.0	H	350.0			
CORRELATION FACTOR		A	1	B	0			
ON-BOARD STABILITY PERIOD		30 DAYS						
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#	AGE L	AGE H
			SEX	YEAR	MONTH	YEAR	MONTH	
			1	#	#	#	#	#
			2	#	#	#	#	#
			3	#	#	#	#	#
			4	#	#	#	#	#
			5	#	#	#	#	#
			6	#	#	#	#	#
			7	#	#	#	#	#
PANIC VALUE		L	#	H	#	UNIT	U/L	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR	FACTOR			
	POINT 1	*		*	OD-L	OD-H		
	POINT 2	*		*				
	POINT 3	*		*				
	POINT 4	*		*				
	POINT 5	*		*				
	POINT 6	*		*				
	POINT 7	*		*				
1-POINT CAL. POINT								
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD				

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TEST NAME	LAC	TYPE	SERUM	OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl
	R2 VOLUME		µl	DILUTION	0	µl
WAVELENGTH	PRI	520		SEC		
METHOD		END				
REACTION SLOP		+				
MEASURING POIN 1	FIRST	0		LAST	27	
MEASURING POIN 2	FIRST			LAST		
LINEARITY		%				
N0-LAG-TIME						
PRE-DILUTION RATE						
	MIN OD			MAX OD		
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50	
	LAST L	-0.100		LAST H	2.50	
DYNAMIC RANGE	L	0.0		H	120.0	
CORRELATION FACTOR	A	1		B	0	
ON-BOARD STABILITY PERIOD		30 DAYS				
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#
			AGE L		AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH	
	1	#	#	#	#	
	2	#	#	#	#	
	3	#	#	#	#	
	4	#	#	#	#	
	5	#	#	#	#	
	6	#	#	#	#	
	7	#	#	#	#	
PANIC VALUE	L	#	H	#	UNIT	mg/dL
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H	
	POINT 1	*		*		
	POINT 2	*		*		
	POINT 3	*		*		
	POINT 4	*		*		
	POINT 5	*		*		
	POINT 6	*		*		
	POINT 7	*		*		
1-POINT CAL. POINT						
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD		

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TEST NAME	LDH		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	4	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD	RATE							
REACTION SLOP		-						
MEASURING POIN 1	FIRST	13		LAST	25			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	2000.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#		H	#	UNIT	U/L	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	LDH		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	4	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD	RATE							
REACTION SLOP		-						
MEASURING POIN 1	FIRST	8		LAST	20			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	2000.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	U/L		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL Nº	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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TEST NAME	LIPASE		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	160	µl	DILUTION	0	µl		
	R2 VOLUME	40	µl	DILUTION	0	µl		
WAVELENGTH	PRI	578		SEC				
METHOD		FIXED						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	9		LAST	17			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	300.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD		30 DAYS						
NORMAL RANGES	VALUE/FLAG	#		LEVE L	#		LEVEL H	#
				AGE L			AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH			
	1	#	#	#	#			
	2	#	#	#	#			
	3	#	#	#	#			
	4	#	#	#	#			
	5	#	#	#	#			
	6	#	#	#	#			
	7	#	#	#	#			
PANIC VALUE	L	#		H	#		UNIT	UI/mL
CALIBRATION TYPE	AB			FORMULA	Y=AX+B		COUNTS	#
	CAL N°	OD	CONC	FACTOR	OD-L	FACTOR	OD-H	
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR				CALIBRATION STABILITY PERIOD				

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TEST NAME	MAG		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	520		SEC				
METHOD		END						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	3.5			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#	H	#	UNIT	mg/dL		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B		COUNTS	#		
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR								
					CALIBRATION STABILITY PERIOD			

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TEST NAME	PCR TURBI		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	1.5	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	240	µl	DILUTION	0	µl		
	R2 VOLUME	60	µl	DILUTION	0	µl		
WAVELENGTH	PRI	540		SEC				
METHOD	FIXED							
REACTION SLOP	+							
MEASURING POIN 1	FIRST	4		LAST	12			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	150.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#		H	#	UNIT	mg/L	
CALIBRATION TYPE	AB		FORMULA	Y=AX+B		COUNTS	#	
	CAL N°	OD	CONC	FACTOR	FACTOR			
POINT 1	*		*	OD-L	OD-H			
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR			CALIBRATION STABILITY PERIOD					

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TEST NAME	PCR ULTRA		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	<input type="text" value="2"/>	µl	DILUTION	<input type="text" value="0"/>	µl		
REAGENTS	R1 VOLUME	<input type="text" value="160"/>	µl	DILUTION	<input type="text" value="0"/>	µl		
	R2 VOLUME	<input type="text" value="40"/>	µl	DILUTION	<input type="text" value="0"/>	µl		
WAVELENGTH	PRI	<input type="text" value="540"/>		SEC	<input type="text"/>			
METHOD		<input type="text" value="FIXED"/>						
REACTION SLOP		<input type="text" value="+"/>						
MEASURING POIN 1	FIRST	<input type="text" value="4"/>		LAST	<input type="text" value="20"/>			
MEASURING POIN 2	FIRST	<input type="text"/>		LAST	<input type="text"/>			
LINEARITY		<input type="text"/>	%					
N0-LAG-TIME		<input type="text"/>						
PRE-DILUTION RATE		<input type="text"/>						
	MIN OD	<input type="text"/>		MAX OD	<input type="text"/>			
REAGENT OD LIMIT	FIRST L	<input type="text" value="-0.100"/>		FIRST H	<input type="text" value="2.50"/>			
	LAST L	<input type="text" value="-0.100"/>		LAST H	<input type="text" value="2.50"/>			
DYNAMIC RANGE	L	<input type="text" value="0.0"/>		H	<input type="text" value="10.0"/>			
CORRELATION FACTOR	A	<input type="text" value="1"/>		B	<input type="text" value="0"/>			
ON-BOARD STABILITY PERIOD		<input type="text" value="30 DAYS"/>						
NORMAL RANGES	VALUE/FLAG	<input type="text" value="#"/>	LEVE L	<input type="text" value="#"/>	LEVEL H	<input type="text" value="#"/>		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
	1	#	#	#	#	#		
	2	#	#	#	#	#		
	3	#	#	#	#	#		
	4	#	#	#	#	#		
	5	#	#	#	#	#		
	6	#	#	#	#	#		
	7	#	#	#	#	#		
PANIC VALUE	L	<input type="text" value="#"/>		H	<input type="text" value="#"/>	UNIT	<input type="text" value="mg/L"/>	
CALIBRATION TYPE		<input type="text" value="AB"/>	FORMULA	<input type="text" value="Y=AX+B"/>	COUNTS	<input type="text" value="#"/>		
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT		<input type="text"/>						
MB TYPE FACTOR		<input type="text"/>		CALIBRATION STABILITY PERIOD	<input type="text"/>			

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TEST NAME	<input type="text" value="PROT"/>	TYPE	<input type="text" value="SERUM"/>	OPERATION	<input type="text" value="YES"/>	
SAMPLE	VOLUME	<input type="text" value="2"/>	µl	DILUTION	<input type="text" value="0"/>	µl
REAGENTS	R1 VOLUME	<input type="text" value="200"/>	µl	DILUTION	<input type="text" value="0"/>	µl
	R2 VOLUME	<input type="text"/>	µl	DILUTION	<input type="text" value="0"/>	µl
WAVELENGTH	PRI	<input type="text" value="540"/>		SEC	<input type="text"/>	
METHOD		<input type="text" value="END"/>				
REACTION SLOP		<input type="text" value="+"/>				
MEASURING POIN 1	FIRST	<input type="text" value="0"/>		LAST	<input type="text" value="27"/>	
MEASURING POIN 2	FIRST	<input type="text"/>		LAST	<input type="text"/>	
LINEARITY		<input type="text"/>	%			
N0-LAG-TIME		<input type="text"/>				
PRE-DILUTION RATE		<input type="text"/>				
	MIN OD	<input type="text"/>		MAX OD	<input type="text"/>	
REAGENT OD LIMIT	FIRST L	<input type="text" value="-0.100"/>		FIRST H	<input type="text" value="2.50"/>	
	LAST L	<input type="text" value="-0.100"/>		LAST H	<input type="text" value="2.50"/>	
DYNAMIC RANGE		L	<input type="text" value="0.0"/>	H	<input type="text" value="12.0"/>	
CORRELATION FACTOR		A	<input type="text" value="1"/>	B	<input type="text" value="0"/>	
ON-BOARD STABILITY PERIOD		<input type="text" value="30 DAYS"/>				
NORMAL RANGES	VALUE/FLAG	<input type="text" value="#"/>	LEVE L	<input type="text" value="#"/>	LEVEL H	<input type="text" value="#"/>
			AGE L		AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH	
	1	#	#	#	#	
	2	#	#	#	#	
	3	#	#	#	#	
	4	#	#	#	#	
	5	#	#	#	#	
	6	#	#	#	#	
	7	#	#	#	#	
PANIC VALUE	L	<input type="text" value="#"/>	H	<input type="text" value="#"/>	UNIT	<input type="text" value="g/dL"/>
CALIBRATION TYPE	<input type="text" value="AB"/>	FORMULA	<input type="text" value="Y=AX+B"/>	COUNTS	<input type="text" value="#"/>	
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H	
POINT 1	*		*			
POINT 2	*		*			
POINT 3	*		*			
POINT 4	*		*			
POINT 5	*		*			
POINT 6	*		*			
POINT 7	*		*			
1-POINT CAL. POINT	<input type="text"/>					
MB TYPE FACTOR	<input type="text"/>			CALIBRATION STABILITY PERIOD	<input type="text"/>	

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TEST NAME	PROTU		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	4	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	200	µl	DILUTION	0	µl		
	R2 VOLUME		µl	DILUTION	0	µl		
WAVELENGTH	PRI	600		SEC				
METHOD		END						
REACTION SLOP		+						
MEASURING POIN 1	FIRST	0		LAST	27			
MEASURING POIN 2	FIRST			LAST				
LINEARITY		%						
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	3000.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#	LEVE L	#	LEVEL H	#		
			AGE L		AGE H			
	SEX	YEAR	MONTH	YEAR	MONTH			
	1	#	#	#	#	#		
	2	#	#	#	#	#		
	3	#	#	#	#	#		
	4	#	#	#	#	#		
	5	#	#	#	#	#		
	6	#	#	#	#	#		
	7	#	#	#	#	#		
PANIC VALUE	L	#	H	#	UNIT	mg/L		
CALIBRATION TYPE	AB	FORMULA	Y=AX+B	COUNTS	#			
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR		CALIBRATION STABILITY PERIOD						

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1-POINT CAL. POINT

MB TYPE FACTOR CALIBRATION STABILITY PERIOD

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TEST NAME	UREIA		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	<input type="text" value="2"/>	µl	DILUTION	<input type="text" value="0"/>	µl		
REAGENTS	R1 VOLUME	<input type="text" value="200"/>	µl	DILUTION	<input type="text" value="0"/>	µl		
	R2 VOLUME	<input type="text"/>	µl	DILUTION	<input type="text" value="0"/>	µl		
WAVELENGTH	PRI	<input type="text" value="340"/>		SEC	<input type="text"/>			
METHOD		<input type="text" value="FIXED"/>						
REACTION SLOP		<input type="text" value="-"/>						
MEASURING POIN 1	FIRST	<input type="text" value="4"/>		LAST	<input type="text" value="12"/>			
MEASURING POIN 2	FIRST	<input type="text"/>		LAST	<input type="text"/>			
LINEARITY		<input type="text"/>	%					
N0-LAG-TIME		<input type="text"/>						
PRE-DILUTION RATE		<input type="text"/>						
	MIN OD	<input type="text"/>		MAX OD	<input type="text"/>			
REAGENT OD LIMIT	FIRST L	<input type="text" value="-0.100"/>		FIRST H	<input type="text" value="2.50"/>			
	LAST L	<input type="text" value="-0.100"/>		LAST H	<input type="text" value="2.50"/>			
DYNAMIC RANGE	L	<input type="text" value="0.0"/>		H	<input type="text" value="300.0"/>			
CORRELATION FACTOR	A	<input type="text" value="1"/>		B	<input type="text" value="0"/>			
ON-BOARD STABILITY PERIOD		<input type="text" value="30 DAYS"/>						
NORMAL RANGES	VALUE/FLAG	<input type="text" value="#"/>		LEVE L	<input type="text" value="#"/>		LEVEL H	<input type="text" value="#"/>
				AGE L			AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH			
	1	#	#	#	#			
	2	#	#	#	#			
	3	#	#	#	#			
	4	#	#	#	#			
	5	#	#	#	#			
	6	#	#	#	#			
	7	#	#	#	#			
PANIC VALUE	L	<input type="text" value="#"/>		H	<input type="text" value="#"/>		UNIT	<input type="text" value="mg/dL"/>
CALIBRATION TYPE		<input type="text" value="AB"/>		FORMULA	<input type="text" value="Y=AX+B"/>		COUNTS	<input type="text" value="#"/>
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
	POINT 1	*		*				
	POINT 2	*		*				
	POINT 3	*		*				
	POINT 4	*		*				
	POINT 5	*		*				
	POINT 6	*		*				
	POINT 7	*		*				
1-POINT CAL. POINT		<input type="text"/>						
MB TYPE FACTOR		<input type="text"/>		CALIBRATION STABILITY PERIOD	<input type="text"/>			

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TEST NAME	UREIA		TYPE	SERUM		OPERATION	YES	
SAMPLE	VOLUME	2	µl	DILUTION	0	µl		
REAGENTS	R1 VOLUME	140	µl	DILUTION	0	µl		
	R2 VOLUME	60	µl	DILUTION	0	µl		
WAVELENGTH	PRI	340		SEC				
METHOD	FIXED							
REACTION SLOP	-							
MEASURING POIN 1	FIRST	10		LAST	18			
MEASURING POIN 2	FIRST			LAST				
LINEARITY	%							
N0-LAG-TIME								
PRE-DILUTION RATE								
	MIN OD			MAX OD				
REAGENT OD LIMIT	FIRST L	-0.100		FIRST H	2.50			
	LAST L	-0.100		LAST H	2.50			
DYNAMIC RANGE	L	0.0		H	300.0			
CORRELATION FACTOR	A	1		B	0			
ON-BOARD STABILITY PERIOD	30 DAYS							
NORMAL RANGES	VALUE/FLAG	#		LEVE L	#		LEVEL H	#
				AGE L			AGE H	
	SEX	YEAR	MONTH	YEAR	MONTH			
1	#	#	#	#	#			
2	#	#	#	#	#			
3	#	#	#	#	#			
4	#	#	#	#	#			
5	#	#	#	#	#			
6	#	#	#	#	#			
7	#	#	#	#	#			
PANIC VALUE	L	#		H	#		UNIT	mg/dL
CALIBRATION TYPE	AB			FORMULA	Y=AX+B		COUNTS	#
	CAL N°	OD	CONC	FACTOR OD-L	FACTOR OD-H			
POINT 1	*		*					
POINT 2	*		*					
POINT 3	*		*					
POINT 4	*		*					
POINT 5	*		*					
POINT 6	*		*					
POINT 7	*		*					
1-POINT CAL. POINT								
MB TYPE FACTOR								
	CALIBRATION STABILITY PERIOD							

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